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These instructions for use must be followed before their clinical use.

Notes regarding the delivered goods.

The individual parts of the system may only be accepted when the manufacturer's label and packaging are undamaged and unopened at the time of delivery. If this is not the case, the rejected goods must be returned to Industrial Médica Alicante S.L., Elche / Spain or its authorized representative within ten working days.

The single-use devices can not be reused, since, intended for use in a single person during a surgical procedure, and are not designed to work as intended after the first use. All components are served **NON-STERILE** and must be prepared properly before their first use. All packaging must be removed prior to being used.

Information about the material

All INNOVA implants are made of grade 2 titanium (ASTM F67, ISO 5832-2) and ELI grade 23 titanium (ASTM F136, ISO 5832-3). All the titanium materials used are biocompatible, anticorrosive, and non-toxic in a biological environment.

Regarding the associated instruments, the Kirschner wires are manufactured from AISI 316vIm stainless steel (ASTM F 138); The rest of the instruments are made of stainless steel, titanium, PPP and aluminium.

Color Code

A color code is used in Implants and instruments, to highlight those that have a specific purpose, and those that are associated with each other:

Gold Screws	Quick Fix, Variable angle locking
Pink Screws	Quick Fix, Variable angle non locking
Blue Screws	SR FIX, cannulated self-drilling compression screw, 2.5mmØ (Scarff osteotomy)

Green Screws	SR FIX, cannulated self-drilling compression screw, 3.5mmØ (Chevron/Austin osteotomy)
Pink Snap-off screws	SR FIX, Snap-off screws 2.0 mmØ
Blue ring instruments	Screwdriver, drill bit and countersink to use with SR FIX 2.5 mmØ screws
Green ring instruments	Screwdriver, drill bit and countersink to use with SR FIX 3.5 mmØ screws
Pink ring instruments	Screwdriver to use with snap-off 2.0 mmØ screws

Intended use

INNOVA F-SYSTEM implants are single-use devices, designed for temporary fixation, correction and / or bone stabilization. The implants have been manufactured from recognized and accepted materials for implantation in the body, thus complying with the ASTM and / or ISO standards:

- Titanium grade 2, for plates, complying with ISO 5832-2 and ASTM F67-13
- Titanium grade 23 ELI for screws, complying with ISO 5832-3 and ASTM F136-13

Contraindications

The studies, the personal training and the professional criteria of the doctor are indispensable for the selection of the device and the indicated treatment.

The conditions that involve a higher risk of failure are, among others:

- Any active infection, or any marked local inflammation in the affected area or the adjacent one.
- Encompassed vascularity that prevents the blood supply in the affected area.
- Alterations in the bone reserve due to a disease, infection or previous implantation, which can not provide adequate support and / or fixation for the devices.
- Sensitivity to the material used.
- Obesity. An overweight or obese patient can produce overloads in the implant that can lead to failures in the fixation of the device or in the device itself.
- Patients with inadequate or deficient tissue coverage.
- Implant interference with anatomical structures or physiological performance.
- Any mental or neuromuscular disorder that generates an unacceptable risk of fixation failure or postoperative complications.

Other medical or surgical conditions that prevent possible surgical benefits.

WARNINGS AND PRECAUTIONS

- Selection and size of the implant:**

IMA, as a manufacturer, does not recommend a specific surgical procedure for each patient. It is very important to select the correct fixation device, so it is the responsibility of the healthcare professional to select the appropriate implant for each clinical case.

In general, implants are designed to remain in the body temporarily and be explanted after sufficient healing (bone) has taken place. The estimated operative period of the implants is between 30 days and 6 months.

The surgeon must guarantee or favor osteosynthesis in terms of adaptation or stability of mobilization (for example, splint and / or immobilization) in the postoperative process, until the bone healing is completed. Early loading and / or restriction of blood

supply may increase the risk of loosening, migration or rupture of the devices. In the case of complications, it may be necessary to remove the implants. For its extraction, use the instruments that I.M.A. supplies. Make sure the screwdriver / screw head connection are aligned.

It will be the patient's decision to retain or have the implants explanted after the treatment is finished.

The adequate blood flow should be maintained and allow a rigid and stable fixation, since it could lead to loosening, bending, cracking or fracture of the device and / or bone. The proper size of the implant for a given patient can be established based on the patient's stature, weight, functional requirements, and anatomy. Each implant must be used in the correct anatomical location, in accordance with accepted internal fixation standards.

Information related to magnetic resonance environments

The safety and compatibility of the devices described in this leaflet in an MR environment have not been evaluated.

Preoperative

The implant is for single use. Make sure that all the components necessary for the intervention are available in the operating room.

It is recommended to perform a pre-surgical inspection to check if the material to be used has suffered damage during storage. Although it is not very common, it is possible that fractures or ruptures of instruments that have been exposed to force and / or excessive use may occur during the surgical act. Therefore, we recommend examining the instruments before the surgical intervention.

Prevent the implants from suffering superficial damage.

- Discard all damaged or improperly handled implants.
- Follow the cleaning and sterilization instructions.
- As far as possible, the implant should be prevented from bending or breaking, as it can reduce its resistance to fatigue and cause breakage when subjected to load.
- During the procedure, repeatedly check the firmness of the connection between the implant and the instruments, or the connection between the instruments, which is necessary for precise positioning and fixation.

Do not use the components of the systems of products of Industrial Médica Alicante S.L., together with the components of systems of other manufacturers, unless otherwise specified (see the manual of surgical techniques).

Intraoperative

These implants are not designed to support the entire weight of the patient or a significant part of the weight for long periods of time. For this reason, warnings and postoperative instructions are of the utmost importance. For this reason, external immobilization (for example, orthopedic appliances or plaster) can be used until bone consolidation is confirmed by radiographs or other procedures.

The implant is a device for prolonged use (between 30 days and 6 months). If the bone consolidation is delayed or does not occur, or the explantation is not carried out, complications may arise, such as rupture or loosening of the implant, or instability of the implant system. It is recommended to periodically perform postoperative examinations (for example, radiographs).

The risk of postoperative complications (for example, the failure of an implant) is higher in obese patients, in patients with mental or neuromuscular disorders. For this reason, an additional individualized postoperative follow-up should be performed.

In addition, After removal of an implant, adequate postoperative control must be performed to prevent the bone from fracturing or fracturing again.

ADVERSE EFFECTS

Generally, the adverse effects are due to clinical reasons and are not related to the implanted device. The following are the most

frequent adverse effects of the use of internal fixation devices:

- These devices may fail due to breakage, deformation and / or loosening, if they are subject to a greater load associated with delayed consolidation or non-consolidation. The internal fixation devices are load distribution devices and are designed to place the fractured bone surface in apposition in order to facilitate its consolidation. If there is a delay in the consolidation or it does not occur, the device may eventually break due to metal fatigue. The loads produced on the device by the load support and the activity level of the patient will determine the longevity of the device. Conditions attributable to non-consolidation, osteoporosis, osteomalacia, diabetes, inhibited revascularization.
- Incorrect alignment between fragments may be the cause of non-consolidation of the bone and / or bending, cracking or even rupture of the device.
- Increased fibrous tissue response around the fracture zone due to unstable comminuted fractures
- Immediate or late infection, both deep and superficial.
- Deep vein thrombosis.
- Avascular necrosis.
- Possible subclinical damage to the nervous system due to surgical trauma.

VA RING® variable angle locking

- Variable angle Correct locking (+ 15°) of the LS screws on the plates with VA RING® hole:**

A visual inspection of the screw head provides an indicator of correct locking.

To ensure optimal locking of the screw on the plate, check that the screw head is fully inserted into the hole in the plate. (Fig. 1). If the protrusion of the screw head is evident (Fig. 2), it has not been fully inserted into the hole in the plate and has not reached the locking position. We recommend the use of the torque limiter Ref. 10.0210 to avoid damaging the plate and ensure locking functionality. In this case, the tightening of the screw must be repeated to obtain a complete insertion and a correct locking. In cases of poor bone quality, a slight axial pressure will help achieve it.

Additional Information

Additional information about the products (for example, the surgical technique, the operating instructions for care, cleaning, disinfection and sterilization of non-sterile products) can be requested from your local consultant at Industrial Médica Alicante S.L. or distributor. In addition, all the relevant information can be found on the Internet at www.industrialmedicaalicantina.com.

HANDLING AND STERILIZATION

Any change in the mechanical, physical or chemical characteristics that has been caused under repeated conditions of use, cleaning and sterilization, may compromise the integrity of the design and / or materials, reducing safety, performance and / or compliance with the relevant specifications. Consult the label of the device to identify if it is one or several uses and / or its removal for cleaning and sterilization

Cleaning

- All implants, instruments and containers that are **NOT STERILE** must be cleaned, disinfected and sterilized before each use. Cleaning and mechanical inspection prior to sterilization must be carried out by qualified personnel.
- Comply with the indications and recommendations regarding chemical detergents in the instructions for use provided by the manufacturers of the cleaning and sterilization equipment (manual and / or machine cleaning, ultrasonic treatment, etc.).

- Do not use an implant if the surface has been damaged. Damaged implants should be discarded.
- All users must be qualified personnel with documentary evidence of training and competence. Users must wear personal protective equipment (PPE).

Cleaning methods

- Manual processing

Equipment: soft brush, neutral enzymatic cleaner or neutral detergent with a 5.5 ≤ pH ≤ 8.5.

- Prepare a solution with warm current water and detergent or cleaner. Follow the manufacturer's recommendations for the detergent or enzymatic cleaner and pay close attention to the correct exposure time, temperature, water quality and concentration.
- Wash the implant carefully and by hand. Do not use steel wool or abrasive cleaners with implants. Rinse the implant thoroughly with deionized or purified water. Use deionized or purified water for final rinsing.
- Dry the implant with a clean, soft, lint-free cloth so as not to scratch the surface.

- Ultrasound processing

Equipment: Ultrasonic cleaner, neutral enzymatic cleaner or neutral detergent with 5.5 ≤ pH ≤ 8.5.

Note: ultrasonic cleaning can cause more damage to implants that already have surface damage.

- Prepare a solution with warm running water and detergent or cleaner. Follow the manufacturer's recommendations for the detergent or enzymatic cleaner and pay close attention to the correct exposure time, temperature, water quality and concentration.
- Clean the implants with ultrasound for a minimum of 15 minutes.
- Thoroughly rinse the implant with deionized or purified water. Use deionized or purified water for final rinsing.
- Dry the implant with a clean, soft, lint-free cloth so as not to scratch the surface.

- Mechanical processing

Equipment: washer / disinfectant, neutral enzymatic cleaner or neutral detergent with a 5.5 ≤ pH ≤ 8.5.

Industrial Médica Alicante S.L. recommends the following instructions for cleaning and sterilizing the instruments:

- Clean all instruments carefully using a mild detergent and brushes and warm water. Ensure that dried blood, bone and other deposits have been completely removed from both the instruments and the sterilization tray.
- Carefully rinse the instruments and the sterilization tray with water.
- Place the instruments in the sterilization box and make sure that the lid is placed and closed correctly.
- Sterilize by autoclaving following the instructions below.

Sterilization

Sterilize by autoclaving with saturated steam, following the instructions below.

- The products are supplied in a **NON-STERILE** condition.
- The use of the sterilization equipment must be done in compliance with the instructions provided by the manufacturer of the sterilization equipment in the instructions for use.
- The products should be cleaned and disinfected at the user's facilities before sterilization following standard hospital procedures.
- Implant re-sterilization should not be performed if the implant comes in contact with any type of contamination (for example, contact with biological tissues, such as body fluids or blood).
- The products are sterilizable by steam (autoclave).

All NON-STERILE products can be sterilized with steam in an autoclave (the equipment must comply with EN 13060 or EN 285 standards).

For initial or consecutive sterilization, Industrial Médica Alicante recommends the following parameters according to the requirements of the usual sterilization standards EN ISO 17665 and ANSI / AAMI ST79:

- Pre-vacuum Steam Sterilization: Temperature MIN.132 ° C and MAX.135 ° C.
- Pre-conditioning pulses: 3 minutes
- Time: ≥ 3-4 min.
- Drying time:> 20 - 30 min.

Configuration: Tray or surgical container with appropriate wrapping.

Packing technique: Double wrap the component in a CE-marked packaging suitable for steam sterilization. Use a sequential technique to wrap the component.

Do not exceed 135 ° C to avoid compromising the functions of the instruments with polymeric components. The sterilization cycle will consist of six main phases according to ISO 17665-1:

- Air extraction phase
- Loading phase
- Sterilization phase.
- Vaporization phase.
- Drying phase
- Equalization phase

Industrial Médica Alicante S.L. recommends sterilization according to the validated methods indicated above. If the user wishes to use other methods (eg flash sterilization), he must validate them. The user will be the final responsible for the validation of the sterilization equipment and techniques.

STORAGE CONDITIONS

Store in a dry, cool place and keep away from direct sunlight. Before use, inspect the product packaging for signs of manipulation or aqueous contamination. Use the oldest lots first.

No special conditions are required for the storage of F-SYSTEM products.

Manufacturer

Industrial Médica Alicante S.L.
Miguel Servet 14
03203 Elche / Spain

	DO NOT REUSE
	USE BY DATE
	LOT NUMBER
	ARTICLE NUMBER
	DATE OF MANUFACTURE
	MANUFACTURER
	CAUTION: CONSULT ACCOMPANYING DOCUMENTS
	SEE INSTRUCTIONS FOR USE
	KEEP AWAY OF SUNLIGHT
	KEEP DRY
	NO ESTÉRIL

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This document is subject to continuous review. Verify that the current printed version is identical to the one in www.industrialmedicaalicantina.com/download/instructions-for-use/.